# K-1 High Shear Lab Mixer

The K-1 Lab Mixer adds a new dimension to the mixing/pelletizing process.

Our high shear lab mixer is ideal for testing, experimentation and development of new processes.

An economical and rugged machine that can also be used for small production work.



### **Specifications**

Attribute	Imperial	Metric
Pan Volume	0.7 cu. ft.	20 L
Working volume	0.3 cu. ft.	8 L
Pan motor	¾ HP	0.6 Kw
Pan speed	Variable to 50 rpm	
Rotor motor	3 HP	2.25 Kw
Rotor speed	Variable to 3000 rpm	
Weight	1100 Lbs	500 Kg
Length	40"	100 cm
Width	30"	75 cm
Height	53"	135 cm
Pan motor Pan speed Rotor motor Rotor speed Weight Length Width	¾ HP  Variable to 50 rpm  3 HP  Variable to 3000 rpm  1100 Lbs  40"  30"	0.6 Kw 2.25 Kv 500 Kg 100 cm 75 cm

## **Design Features:**

- Removable stainless-steel mixing pan with 0.3 cu. ft. capacity
- Large viewing port on pan cover for adding process materials, observing the mix, and sampling
- Touch screen control on independent control stand
- Thermocouple mix temperature reading
- Dual hand safety buttons for operating mixing head

OPTIONS: pan heating and cooling, various rotor styles

## The K-1 Lab Mixer has the same counter-current mixing design as all Lancaster Mixers which provides improved mix quality and can pelletize in one step.




### **HOW THE LANCASTER K1-LAB MIXER WORKS**

The K1-Lab Mixer has the same counter-current mixing design as all Lancaster Mixers with the clockwise rotating pan bringing material into the counterclockwise high speed rotor.

This mixing design provides not only improved mix quality but can pelletize in one process step. The rotor and pan speed are infinitely variable.

These features, combined with the rugged construction and the various styles of rotors available, make the K1-Lab the mixer of choice for almost any type of batch mixing.





#### **ELECTRICAL CONTROL FEATURES**

- Electronic individually controlled speeds of pan and rotor
- Pre-programmable mix design times
- Input or adjust pan and rotor speeds before or during operation
- Data record of operating functions can be exported to a USB or other device



920 Mechanic Street, Lebanon, PA 17046 800.447.7351 www.LancasterProducts.com